

**Recipe Suggestion**  
**Palsgaard® DATEM 3519**  
**Palsgaard® 6111**

**Liquid shortening for yeast raised products**



**Product:** Liquid Shortening for yeast raised products

**Recipe:**

	%
Palsgaard® DATEM 3519	8,00
Palsgaard® 6111	2,00
Liquid rape seed oil	90,0
	100.00

**Procedure:**

- \* Melt the oil/fat blend and adjust its temperature to approx. 45°C
- \* Dissolve Palsgaard® DATEM 3519 and Palsgaard® 6111 in 5 parts of oil/fat at approx. 70°C.
- \* Transfer the emulsifier/oil solution into the oil/fat mix.
- \* Emulsification temperature: approx. 45°C.

Identification: Liquid Shortening for yeast raised products - trans free (00351901-EU-E-RS).doc